

# **Recombinant green fluorescent protein rTurboGFP**

**Cat.# FP552**


<b>Amount:</b>	<b>100 ug</b>
<b>Concentration:</b>	<b>1.0 mg/ml</b>
<b>Storage buffer:</b>	<b>1x PBS</b>
<b>Storage conditions:</b>	<b>store at +4°C</b>
<b>Shipping conditions:</b>	<b>ambient temperature</b>
<b>Shelf life:</b>	<b>6 months from date of receipt under proper storage conditions</b>
<b>Lot number:</b>	<b>specified on product label</b>

Recombinant TurboGFP (rTurboGFP) is 27-kDa green fluorescent protein. It has spectral properties identical to those of the expressed TurboGFP and is suitable as control reagent for TurboGFP expression using the TurboGFP expression vectors.

rTurboGFP is purified from transformed *E. coli* using organic extraction and hydrophobic chromatography or metal-ion affinity chromatography (methods vary for different lots). Both methods ensure high purity of the recombinant protein and maintenance of fluorescence. The protein concentration is measured by chromophore absorption. rTurboGFP contains 6xHis tag at its N-terminus.

**TurboGFP normalized excitation (thin line) and emission (thick line) spectra.**

Complete TurboGFP spectra in Excel format can be downloaded from the Evrogen Web site at [www.evrogen.com/TurboGFP.shtml](http://www.evrogen.com/TurboGFP.shtml)



### **rTurboGFP as standard on protein gel**

As a standard on protein gel, the recombinant protein can be used to correlate TurboGFP expression levels to fluorescence intensity or to differentiate problems with detection of TurboGFP fluorescence from expression of TurboGFP protein.

When denatured by heating (2-3 min, 95-98°C), rTurboGFP is detected as a doublet at about 27 kDa on Coomassie stained SDS gel. If rTurboGFP is to be used as an internal standard in a Coomassie-stained minigel, we recommend loading 0.5-1.0 µg of rTurboGFP per lane. If rTurboGFP is added to a total cell/tissue lysate or other crude sample, the amount of total protein loaded per lane must be optimized for the particular application.

Predenatured rTurboGFP will not fluoresce on SDS gel, whereas nondenatured one will fluoresce on SDS gel.

### **rTurboGFP as control for fluorescence microscopy**

The following protocols are for rTurboGFP use as a control on microscope slides in fluorescence microscopy. The purified proteins may be used to optimize lamp and filter set conditions

for detection of TurboGFP fluorescence, or as a qualitative means to correlate TurboGFP fluorescence with protein amount in transfected cells.

### **A. Unfixed samples**

Please use this method for live cell fluorescence or other cases where a fixation step is not desired

1. Perform 1:10 serial dilutions of the 1.0 mg/ml rTurboGFP stock solution with 10 mM Tris-HCl (pH 8.0) to yield concentrations of 0.1 mg/ml and 0.01 mg/ml.

#### **Notes:**

- These dilutions should suffice as a positive control. The 1.0 mg/ml solution will give a very bright fluorescent signal by microscopy.
  - The diluted samples can be stored at +4°C for up to 3 months with no loss of fluorescence intensity.
2. Using a micropipette, spot 1-2  $\mu$ l of diluted protein onto the microscope slide. If slide contains a mounted coverslip, position the spot several millimeters away from the sample such that a second coverslip can be added over the protein spot.
3. Allow the protein to air-dry for a few seconds, and mark the position of the spot on the other side of the slide to aid in focusing.
4. Add a coverslip over the spot using a 90% glycerol solution in 100 mM Tris-HCl (pH 7.5).
5. Fluorescence from the spot is best viewed at low magnification, using either a 10X or 20X objective lens.

### **B. Fixed samples**

In some cases it may be necessary to fix the recombinant protein to the microscope slide prior to microscopy. This can be

done by dipping the section of the microscope slide containing the air-dried protein spot (after Step A3 above) into 100% methanol for 1 min. Allow the slide to dry completely and place a coverslip over the sample as in Step A4 above.

---

#### **Notice to Purchaser:**

Evrogen Fluorescent Protein Products (the Products) are intended for research use only. The Products are covered by U.S. Pat. # 7,417,131 and other Evrogen Patents and/or Patent applications pending. By use of these Products, you accept the terms and conditions of the applicable Limited Use Label License.

**MATERIAL SAFETY DATA SHEET INFORMATION:** MSDS information is available at: [www.evrogen.com/MSDS.shtml](http://www.evrogen.com/MSDS.shtml)