





Restriction enzyme sites: BseY1, SfiI, BglI, NcoI, MscI

1001 TCCCATCCAACGGCCCTGTGATGCAGAAGAAAACACTCGGCTGGGAGGCCAGCACCCGAGATGCTGTACCCCGCTGACAGCGGCCTGAGAGGCCATGGCCA 1100  
AGGGTAGGTTGCCGGGACACTACGTCTTCTTTTGTGAGCCGACCCTCCGGTCGTGGCTCTACGACATGGGGCGACTGTGCGCCGGACTCTCCGGTACCGGT  
**TurboFP602 > P S N G P V M Q K K T L G W E A S T E M L Y P A D S G L R G H G Q**

Restriction enzyme sites: BssS1, BsgI, Eco57, BpuE1, Bts1, BfuA1, BspM1

1101 GATGGCCCTGAAGCTCGTGGGCGGGGGCTACCTGCACTGCTCCCTCAAGACCACATACAGATCCAAGAAACCCGCTAAGAACCTCAAGATGCCCGGCTTC 1200  
CTACCGGACTTCGAGCACCCGCCCCGATGGACGTGACGAGGGAGTTCTGGTGTATGTCTAGGTTCTTTGGGCGATTCTTGAGTTCTACGGGCCGAAG  
**TurboFP602 > M A L K L V G G G Y L H C S L K T T Y R S K K P A K N L K M P G F**

Restriction enzyme sites: PflM1, Dra3, SpAcc, PflM1, BsaI, PshA1, BssS1, MscI, Scal

1201 CACTTCGTGGACCACAGACTGGAAAAGAATCAAGGAGGCCGACAAAGAGACCTACGTCGAGCAGCACGAGATGGCTGTGGCCAAGTACTGCGACCTCCCTA 1300  
GTGAAGCACCTGGTGTCTGACCTTTCTTAGTTCTCCTCCGGCTGTTTCTCTGGATGCAGCTCGTTCGTGCTCTACCGACACCGGTTTCATGACGCTGGAGGGAT  
**TurboFP602 > H F V D H R L E R I K E A D K E T Y V E Q H E M A V A K Y C D L P S**


Restriction enzyme sites: SpAcc, Bmr1, BspE1, BpuE1, Bgl12, XhoI, Hind3, Sac1, BstB1, EcoR1, SalI, Pst1, Acc65, Kpn1, PspOM, Sac2, BamH1, SmaI, ApaI, XbaI

1301 GCAAAC TGGGGCACAGCTCCGGACTCAGATCTCGAGCTCAAGCTTTCGAATTCTGCAGTCGACGGTACCGCGGGCCCGGGATCCACCGGATCTAGATAACT 1400  
CGTTTGACCCCGTGTCCGAGGCCTGAGTCTAGAGCTCGAGTTTCGAAGCTTAAGACGTCAGCTGCCATGGCGCCCCGGGCCCTAGGTGGCCTAGATCTATTGA  
**TurboFP602/MCS K L G H S S G L R S R A Q A S N S A V D G T A G P G S T G S R \***















```
          NsiI  
        BfrB1 |  
          | |  
4701 CCGTATTACCGCCATGCAT 4719  
      |||||*|||||  
      GGCATAATGGCGGTACGTA
```

**Found:**

Aat2	<b>Acc65</b>	<b>Afe1</b>	<b>Age1</b>	<b>Ale1</b>	AlwN1	<b>Apa1</b>	ApaL1	<b>Ase1</b>	Avr2	<b>BamH1</b>	<b>Bbs1</b>	BciV1	Bcl1
BfrB1	BfuA1	Bgl1	<b>Bgl2</b>	Bmr1	Bpm1	BpuE1	Bsa1	BsaB1	<b>BsaXa</b>	<b>BsaXb</b>	<b>BseR1</b>	BseY1	<b>Bsg1</b>
Bsm1	<b>BspE1</b>	BspH1	<b>BspLU</b>	BspM1	BsrB1	<b>BsrD1</b>	<b>BsrG1</b>	BssS1	BstAP	BstB1	<b>BstX1</b>	Bsu36	BtgZ1
Bts1	<b>Clal</b>	Dra1	Dra3	Drd1	<b>Eag1</b>	Ear1	Eci1	Eco57	<b>EcoR1</b>	<b>Fsp1</b>	<b>Hind3</b>	<b>Hpa1</b>	<b>Kas1</b>
<b>Kpn1</b>	<b>Mfe1</b>	<b>Mlu1</b>	Msc1	Nae1	<b>Nar1</b>	Nco1	<b>Nde1</b>	NgoM4	<b>Nhe1</b>	Nsi1	<b>PflF1</b>	PflM1	polyA
<b>PshA1</b>	Psi1	<b>PspOM</b>	<b>Pst1</b>	Pvu2	<b>Rsr2</b>	<b>Sac1</b>	<b>Sac2</b>	<b>Sall1</b>	Sap1	<b>Sca1</b>	<b>SexA1</b>	Sfi1	<b>Sma1</b>
<b>SnaB1</b>	SpAcc	SpDon	Sph1	Ssp1	<b>Stu1</b>	<b>Xba1</b>	<b>Xho1</b>						

**Unique:**

<b>Acc65</b>	<b>Afe1</b>	<b>Age1</b>	<b>Ale1</b>	<b>Apa1</b>	<b>Ase1</b>	<b>BamH1</b>	<b>Bbs1</b>	<b>Bgl2</b>	<b>BsaXa</b>	<b>BsaXb</b>	<b>BseR1</b>	<b>Bsg1</b>	<b>BspE1</b>
<b>BspLU</b>	<b>BsrD1</b>	<b>BsrG1</b>	<b>BstX1</b>	<b>Clal</b>	<b>Eag1</b>	<b>EcoR1</b>	<b>Fsp1</b>	<b>Hind3</b>	<b>Hpa1</b>	<b>Kas1</b>	<b>Kpn1</b>	<b>Mfe1</b>	<b>Mlu1</b>
<b>Nar1</b>	<b>Nde1</b>	<b>Nhe1</b>	<b>PflF1</b>	<b>PshA1</b>	<b>PspOM</b>	<b>Pst1</b>	<b>Rsr2</b>	<b>Sac1</b>	<b>Sac2</b>	<b>Sall1</b>	<b>Sca1</b>	<b>SexA1</b>	<b>Sma1</b>
<b>SnaB1</b>	<b>Stu1</b>	<b>Xba1</b>	<b>Xho1</b>										

**Not found:**

Aar1	Ac11	Afl2	Ahd1	Asc1	AsiS1	Baela	Baelb	BbvC1	Bcg1a	Bcg1b	Blp1	BmgB1	Bpu10
BsiW1	BsmB1	BssH2	BstE2	BstZ1	_Chi	EcoK	EcoN1	EcoRV	ScFRT	Fse1	FspA1	I_Ceu	loxP
Not1	Nru1	Pac1	Pme1	Pml1	Pvu1	SanD1	Sbf1	Sgf1	SgrA1	Spe1	Srf1	Swal	T3RNA
T7RNA	T7Ter	PISce	Xcm1	Xmn1									

**Excluded by site complexity:**

Acc1	Ac11	Afl3	Alu1	Alw1	Apo1	Ava1	Ava2	Ban1	Ban2	Bbv1	BceA1	Bfal	Bme15
BsaA1	BsaH1	BsaJ1	BsaW1	BseM2	BsiE1	BsiH1	Bsl1	BsmA1	BsmF1	Bsp12	BspCa	BspCb	Bsr1
BsrF1	BssK1	BstF5	BstN1	BstU1	BstY1	Btg1	Cac8	CviJ1	Dde1	Eae1	EcoO1	Fau1	Fnu4H
Fok1	Hae2	Hae3	Hga1	Hha1	Hinc2	Hinf1	HinP1	Hpa2	Hph1	Hpy99	Hpy1	Hpy3	HpyC3
HpyC4	HpyC5	Mae3	Mbo2	Mnl1	Mse1	Msl1	MspA1	Mwo1	Nci1	Nla3	Nla4	Nsp1	Ple1
PpuM1	Rsa1	Sau3A	Sau96	SfaN1	Sfc1	Sml1	Sty1	Taq1	Tat1	Tfi1	Tse1	Tsp45	Tsp50
TspR1													