

pFusionRed-N vector restriction map

The data has not been verified by restriction digestion with each enzyme listed and does not take into account possible methylation sites. Enzymes that recognize unambiguous sequences less than 6 basepairs long are not included – for the more complete enzyme list please refer to the Table of restriction sites.

Unique sites are shown in bold blue. The location given specifies the 3' end of the cut DNA (the base to the left of the cut site). MCS sequence is shown in frame, amino acids encoded by MCS are shown in black.




```

                                     Bsu36      SpDon      Bbs1      Bpm1
                                     |          |          |          |
CCGGGCATCCCCGACTTCTTTAAGCAGTCCTTCCCTGAGGGCTTCACATGGGAGAGAGTCACCACATACGAAGACGGGGCGTGCTGACCGCTACCCAGG
901  |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 1000
GGCCCGTAGGGGCTGAAGAAATTCGTCAGGAAGGGACTCCCGAAGTGTACCCTCTCTCAGTGGTGTATGCTTCTGCCCCCGCACGACTGGCGATGGGTCC
FusionRed > P G I P D F F K Q S F P E G F T W E R V T T Y E D G G V L T A T Q D

                                     BseY1      BseY1
                                     |          |
ACACCAGCCTCCAGGACGGCTGCCTCATCTACAACGTCAAGGTTAGAGGGGTGAACTTCCCAGCCAACGGCCCTGTGATGCAGAAGAAAACACTCGGCTG
1001 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 1100
TGTGGTTCGGAGGTCTGCGGACGGAGTAGATGTTGCAGTTCCAATCTCCCCACTTGAAGGGTTCGGTTGCCGGGACACTACGTCTTCTTTTGTGAGCCGAC
FusionRed > T S L Q D G C L I Y N V K V R G V N F P A N G P V M Q K K T L G W

      BsmB1
StuI  |  Bgl1      BssS1      Eco57      Bsa1
|      |          |          |          |
GGAGGCCTCCACCGAGACGATGTACCCCGCTGACGGCGGCCTGGAAGGCGCATGTGACATGGCCCTGAAGCTCGTGGGCGGGGGCCACCTGATCTGCAAC
1101 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 1200
CCTCCGGAGGTGGCTCTGCTACATGGGGCGACTGCCGCCGGACCTTCCGCGTACACTGTACCGGGACTTCGAGCACCCGCCCGGTGGACTAGACGTTG
FusionRed > E A S T E T M Y P A D G G L E G A C D M A L K L V G G G H L I C N

      BpuE1      SpAcc      Bsa1
      BpuE1      |          PflM1      |
      |          |          |          |
CTTGAGACCACATACAGATCCAAGAAACCCGCTACGAACCTCAAGATGCCCGGCGTCTACAACGTGGACCACAGACTGGAAAGAATCAAGGAGGCCGACG
1201 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 1300
GAACTCTGGTGTATGTCTAGGTTCTTTGGGCGATGCTTGGAGTTCTACGGGCCGAGATGTTGCACCTGGTGTCTGACCTTTCTTAGTTCCCTCCGGCTGC
FusionRed > L E T T Y R S K K P A T N L K M P G V Y N V D H R L E R I K E A D D

                                     Not1
                                     Eag1
      PshA1      BssS1 Dra3      Msc1      SpDon      BsrB1      Xba1      BsaB1
      |          |          |          |          |          |          |
ATGAGACCTACGTCGAGCAGCACGAGGTGGCTGTGGCCAGATACTCTACTGGTGGCGCTGGTGTATGGAGGTAATGAGCGGCCGCGACTCTAGATCATAA
1301 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 1400
TACTCTGGATGCAGCTCGTCGTGCTCCACCGACACCGGTCTATGAGATGACCACCGCGACCACTACCTCCATTTACTCGCCGGCGCTGAGATCTAGTATT
FusionRed > E T Y V E Q H E V A V A R Y S T G G A G D G G K *
```



```

                                polyA      polyA
                                |           |
GGAGACAATACCGGAAGGAACCCGCGCTATGACGGCAATAAAAAGACAGAATAAAACGCACGGTGTGGGTGCGTTTGTTCATAAACGCGGGGTTTCGGTCC
3601 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 3700
CCTCTGTTATGGCCTTCCTTGGGCGCGATACTGCCGTTATTTTTCTGTCTTATTTTGCCTGCCACAACCCAGCAAACAAGTATTTGCGCCCCAAGCCAGG

                                Bsa1
                                |
CAGGGCTGGCACTCTGTTCGATACCCACCGAGACCCCATTTGGGGCCAATACGCCCGGTTTCTTCTTTTTCCCAACCCACCCCAAGTTTCGGGTGAAG
3701 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 3800
GTCCCGACCGTGAGACAGCTATGGGGTGGCTCTGGGGTAACCCCGGTTATGCGGGCGCAAAGAAGGAAAAGGGGTGGGGTGGGGGTTCAAGCCCACTTC

                                BstAP      AlwN1      Bsu36      DraI      DraI
                                |           |           |           |           |
GCCCAGGGCTCGCAGCCAACGTCGGGGCGGCAGGCCCTGCCATAGCCTCAGGTTACTCATATATACTTTAGATTGATTTAAAACCTTCATTTTTAATTTAA
3801 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 3900
CGGGTCCCGAGCGTCGGTTGCAGCCCCGCCGTCCGGGACGGTATCGGAGTCCAATGAGTATATATGAAATCTAACTAAATTTTGAAGTAAAAATTAATTT

                                BspH1
                                |
AAGGATCTAGGTGAAGATCCTTTTTGATAATCTCATGACCAAAATCCCTTAACGTGAGTTTTCGTTCCACTGAGCGTCAGACCCCGTAGAAAAGATCAAA
3901 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 4000
TTCCTAGATCCACTTCTAGGAAAAACTATTAGAGTACTGGTTTTAGGGAATTGCACTCAAAGCAAGGTGACTCGCAGTCTGGGGCATCTTTTCTAGTTT

                                BpuE1
                                |
GGATCTTCTTGAGATCCTTTTTTTCTGCGCGTAATCTGCTGCTTGCAAACAAAAAAACCACCGCTACCAGCGGTGGTTTGTGTTGCCGGATCAAGAGCTAC
4001 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 4100
CCTAGAAGAACTCTAGGAAAAAAGACGCGCATTAGACGACGAACGTTTGTTTTTTTGGTGGCGATGGTGCACCAAAACAAACGGCCTAGTTCTCGATG

                                Eco57      SpAcc
                                |           |
CAACTCTTTTCCGAAGGTAACCTGGCTTCAGCAGAGCGCAGATACCAAATACTGTCTTCTAGTGTAGCCGTAGTTAGGCCACCACTTCAAGAACTCTGT
4101 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 4200
GTTGAGAAAAAGGCTTCCATTGACCGAAGTCGTCTCGCGTCTATGGTTTATGACAGGAAGATCACATCGGCATCAATCCGGTGGTGAAGTTCTTGAGACA

```



```

                                     AlwN1                               BpuE1
AGCACCGCCTACATACCTCGCTCTGCTAATCCTGTTACCAGTGGCTGCTGCCAGTGGCGATAAGTCGTGTCTTACCGGGTTGGACTCAAGACGATAGTTA
4201 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 4300
TCGTGGCGGATGTATGGAGCGAGACGATTAGGACAATGGTCACCGACGACGGTCACCGCTATTTCAGCACAGAATGGCCCAACCTGAGTTCTGCTATCAAT

                                     ApaL1       BseY1                                               SpAcc
CCGGATAAGGCGCAGCGGTTCGGGCTGAACGGGGGGTTCGTGCACACAGCCCAGCTTGGAGCGAACGACCTACACCGAACTGAGATACCTACAGCGTGAGC
4301 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 4400
GGCCTATTCCGCGTCGCCAGCCGACTTGCCCCCAAGCACGTGTGTGGGTTCGAACCTCGCTTGCTGGATGTGGCTTGACTCTATGGATGTGCACTCG

                                     BciV1
                                     EciI       BssS1
TATGAGAAAGCGCCACGCTTCCCGAAGGGAGAAAGGCGGACAGGTATCCGGTAAGCGGCAGGGTCGGAACAGGAGAGCGCACGAGGGGAGCTTCCAGGGGG
4401 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 4500
ATACTCTTTCGCGGTGCGAAGGGCTTCCCTCTTTCGCGCTGTCCATAGGCCATTCGCCGTCCCAGCCTTGTCTCTCGCGTGCTCCCTCGAAGGTCCCC

                                     SpAcc       Drd1       BpuE1       SpAcc       EciI
AAACGCCTGGTATCTTTATAGTCCTGTTCGGTTTTCGCCACCTCTGACTTGAGCGTCGATTTTTGTGATGCTCGTCAGGGGGGCGGAGCCTATGGAAAAAC
4501 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 4600
TTTGCGGACCATAGAAATATCAGGACAGCCAAAGCGGTGGAGACTGAACTCGCAGCTAAAAAACTACGAGCAGTCCCCCGCCTCGGATACCTTTTTTG

                                     SpDon
                                     BspLU
GCCAGCAACCGGCCTTTTTACGGTTCCTGGCCTTTTGCTGGCCTTTTGCTCACATGTTCTTTCCTGCGTTATCCCCTGATTCTGTGGATAACCGTATTA
4601 |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||* 4700
CGGTCGTTGCGCCGAAAAATGCCAAGGACCGGAAAAACGACCGGAAAAACGAGTGTACAAGAAAGGACGCAATAGGGGACTAAGACACCTATTGGCATAAT

                                     NsiI
                                     BfrB1
CCGCCATGCAT
4701 |||||*| 4711
GGCGGTACGTA
```

Found:

Aat2	Acc65	Afe1	Afl2	Age1	Ale1	AlwN1	Apa1	ApaL1	Ase1	Avr2	BamH1	Bbs1	BciV1
BfrB1	Bgl1	Bgl2	Bmr1	Bpm1	BpuE1	Bsa1	BsaB1	BsaXa	BsaXb	BseR1	BseY1	Bsm1	BsmB1
BspH1	BspLU	BspM1	BsrB1	BsrD1	BsrG1	BssS1	BstAP	BstB1	Bsu36	BtgZ1	Bts1	Clal	Dra1
Dra3	Drd1	Eag1	Ear1	Ecil	Ecl2	Eco57	EcoR1	Fsp1	Hind3	Hpa1	Kas1	Kpn1	Mfe1
Msc1	Nae1	Nar1	Nco1	Nde1	NgoM4	Nhe1	Not1	Nsil	PflF1	PflM1	polyA	PshA1	Psi1
PspOM	Pst1	Pvu2	Rsr2	Sac1	Sac2	Sall	Sap1	SexA1	Sfi1	Sma1	SnaB1	SpAcc	SpDon
Sph1	Ssp1	Stu1	Xba1	Xho1									

Unique:

Acc65	Afe1	Afl2	Age1	Ale1	Apa1	Ase1	BamH1	Bbs1	Bgl2	BsaXa	BsaXb	BseR1	BsmB1
BspLU	BsrD1	BsrG1	Bts1	Clal	Ecl2	EcoR1	Fsp1	Hind3	Hpa1	Kas1	Kpn1	Mfe1	Nar1
Nde1	Nhe1	Not1	PflF1	PflM1	PshA1	PspOM	Pst1	Rsr2	Sac1	Sac2	Sall	SexA1	Sfi1
Sma1	SnaB1	Xba1	Xho1										

Not found:

Aar1	Ac11	Ahd1	Asc1	AsiS1	Baela	Baelb	BbvC1	BcglA	Bcglb	Bcl1	Blp1	BmgB1	Bpu10
Bsg1	BsiW1	BspE1	BssH2	BstE2	BstX1	BstZ1	BxatB	BxatL	BxatR	BxatP	_Chi	EcoK	EcoN1
EcoRV	FCatB	FCatL	FCatR	FCatP	ScFRT	Fse1	FspA1	I_Ceu	loxP	Mlu1	Nru1	Pac1	Pme1
Pml1	Pvu1	R4atB	R4atL	R4atP	R4atR	SanD1	Sbf1	Sca1	Sgf1	SgrA1	Spe1	Srf1	Swal
T3RNA	T7RNA	T7Ter	PISce	Xcm1	Xmn1								

Excluded by site complexity:

Acc1	Ac11	Afl3	Alu1	Alw1	Apo1	Ava1	Ava2	Ban1	Ban2	Bbv1	BceA1	Bfa1	Bme15
BsaA1	BsaH1	BsaJ1	BsaW1	BseM2	BsiE1	BsiH1	Bsl1	BsmA1	BsmF1	Bsp12	BspCa	BspCb	Bsr1
BsrF1	BssK1	BstF5	BstN1	BstU1	BstY1	Btg1	Cac8	CviJ1	Dde1	Eae1	EcoO1	Fau1	Fnu4H
Fok1	Hae2	Hae3	Hga1	Hha1	Hinc2	Hinf1	HinP1	Hpa2	Hph1	Hpy99	Hpy1	Hpy3	HpyC3
HpyC4	HpyC5	Mae3	Mbo2	Mnl1	Mse1	Msl1	MspA1	Mwo1	Nci1	Nla3	Nla4	Nsp1	Ple1
PpuM1	Rsa1	Sau3A	Sau96	SfaN1	Sfc1	Sml1	Sty1	Taq1	Tat1	Tfil	Tse1	Tsp45	Tsp50
TspR1													