


```

                Dra1           Dra1           BspH1
                |             |             |
4801 TACTCATATATACTTTAGATTGATTTAAAAC TTCATTTTTAATTTAAAAGGATCTAGGTGAAGATCCTTTTTGATAATCTCATGACCAAAATCCCTTAAC 4900
      |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*
      ATGAGTATATATGAAATCTAACTAAATTTTGAAGTAAAAATTAATTTTCCTAGATCCACTTCTAGGAAAAACTATTAGAGTACTGGTTTTAGGGAATTG

                BpuE1
                |
4901 GTGAGTTTTTCGTTCCACTGAGCGTCAGACCCCGTAGAAAAGATCAAAGGATCTTCTTGAGATCCTTTTTTCTGCGCGTAATCTGCTGCTTGCAAACAAA 5000
      |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*
      CACTCAAAAGCAAGGTGACTCGCAGTCTGGGCATCTTTTCTAGTTTCTAGAAAGAACTCTAGGAAAAAAGACGCGCATTAGACGACGAACGTTTGT

                Eco57
                |
5001 AAAACCACCGCTACCAGCGGTGGTTTGTGTTGCCGGATCAAGAGCTACCAACTCTTTTTCCGAAGGTAAGTGGCTTCAGCAGAGCGCAGATACCAAATACT 5100
      |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*
      TTTTGGTGGCGATGGTCGCCACCAAACAAACGGCCTAGTTCTCGATGGTTGAGAAAAAGGCTTCCATTGACCGAAGTCGTCTCGCGTCTATGGTTTATGA

                SpAcc           AlwN1
                |             |
5101 GTCCTTCTAGTGTAGCCGTAGTTAGGCCACCACTTCAAGAACTCTGTAGCACC GCCTACATACCTCGCTCTGCTAATCCTGTTACCAGTGGCTGCTGCCA 5200
      |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*
      CAGGAAGATCACATCGGCATCAATCCGGTGGTGAAGTTCTTGAGACATCGTGCCGGATGTATGGAGCGAGACGATTAGGACAATGGTCACCGACGACGGT

                BpuE1           ApaL1           BseY1
                |             |             |
5201 GTGGCGATAAGTCGTGTCTTACCGGGTTGGACTCAAGACGATAGTTACCGGATAAGGCGCAGCGGTTCGGGCTGAACGGGGGGTTCGTGCACACAGCCCAG 5300
      |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*
      CACCGCTATTCAGCACAGAATGGCCCAACCTGAGTTCTGCTATCAATGGCCTATTCCGCGTCGCCAGCCCAGACTTGCCCCCAAGCACGTGTGTGCGGGTC

                SpAcc           Eci1
                |             |
5301 CTTGGAGCGAACGACCTACACCGAACTGAGATACCTACAGCGTGAGCTATGAGAAAGCGCCACGCTTCCCGAAGGGAGAAAGCGGACAGGTATCCGGTA 5400
      |||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*|||||*
      GAACCTCGCTTGCTGGATGTGGCTTGACTCTATGGATGTCGCACTCGATACTCTTTCGCGGTGCGAAGGGCTTCCCTCTTTCCGCTGTCCATAGGCCAT
```


Found:

Aar1	Aat2	Afe1	Age1	Ahd1	Ale1	AlwN1	ApaL1	Ase1	Avr2	BamH1	Bbs1	BciV1	Bcl1
BfrB1	BfuA1	Bgl1	Bgl2	BmgB1	Bmr1	Bpm1	BpuE1	Bsa1	BsaB1	BsaXa	BsaXb	BseR1	BseY1
Bsm1	BspE1	BspH1	BspLU	BspM1	BsrB1	BsrD1	BsrG1	BssS1	BstAP	BstB1	Bsu36	BtgZ1	Bts1
_Chi	Clal	Dra1	Dra3	Drd1	Eag1	Ear1	Eci1	Ecl2	Eco57	EcoN1	EcoR1	EcoRV	Fsp1
Hind3	Hpa1	Kas1	Mfe1	Mlu1	Msc1	Nae1	Nar1	Nco1	Nde1	NgoM4	Nhe1	Nsi1	PflF1
polyA	PshA1	Psi1	Pvu2	Rsr2	Sac1	Sap1	Scal	SexA1	Sfil	SnaB1	SpAcc	SpDon	Sph1
Ssp1	Stu1	Xba1	Xcm1	Xho1									

Unique:

Aar1	Afe1	Age1	Ahd1	Ale1	Ase1	BamH1	Bcl1	Bgl2	BmgB1	BspE1	BsrD1	BsrG1	Bts1
_Chi	Eag1	Ecl2	EcoR1	Fsp1	Hind3	Hpa1	Kas1	Mfe1	Mlu1	Nar1	Nde1	Nhe1	PflF1
PshA1	Rsr2	Sac1	Scal	SexA1	Sfil	SnaB1	Xba1	Xcm1					

Not found:

Acc65	Acl1	Afl2	Apa1	Asc1	AsiS1	Bae1a	Bae1b	BbvC1	Bcg1a	Bcg1b	Blp1	Bpu10	Bsg1
BsiW1	BsmB1	BssH2	BstE2	BstX1	BstZ1	BxatB	BxatL	BxatR	BxatP	EcoK	FCatB	FCatL	FCatR
FCatP	ScFRT	Fse1	FspA1	I_Ceu	Kpn1	loxP	Not1	Nru1	Pac1	PflM1	Pme1	Pml1	PspOM
Pst1	Pvu1	R4atB	R4atL	R4atP	R4atR	Sac2	Sall	SanD1	Sbf1	Sgf1	SgrA1	Sma1	Spe1
Srf1	Swa1	T3RNA	T7RNA	T7Ter	PISce	Xmn1							

Excluded by site complexity:

Acc1	Acil	Afl3	Alu1	Alw1	Apo1	Ava1	Ava2	Ban1	Ban2	Bbv1	BceA1	Bfa1	Bme15
BsaA1	BsaH1	BsaJ1	BsaW1	BseM2	BsiE1	BsiH1	Bs11	BsmA1	BsmF1	Bsp12	BspCa	BspCb	Bsr1
BsrF1	BssK1	BstF5	BstN1	BstU1	BstY1	Btg1	Cac8	CviJ1	Dde1	Eae1	EcoO1	Fau1	Fnu4H
Fok1	Hae2	Hae3	Hga1	Hha1	Hinc2	Hinf1	HinP1	Hpa2	Hph1	Hpy99	Hpy1	Hpy3	HpyC3
HpyC4	HpyC5	Mae3	Mbo2	Mnl1	Mse1	Msl1	MspA1	Mwo1	Nci1	Nla3	Nla4	Nsp1	Ple1
PpuM1	Rsa1	Sau3A	Sau96	SfaN1	Sfc1	Sml1	Sty1	Taq1	Tat1	Tfi1	Tse1	Tsp45	Tsp50
TspR1													